## Gabriel

### Gabriel internal test report for bleach cleanability

Test performed:	12 August 2024
Test:	BIFMA HCF 8.1-2019 Health Care Furniture design guidelines or cleanability & ACT Test Method 1-2020
Bleach concentration:	1:10 Sodium Hypochlorite 5.25 – 6.25 %
Product tested:	2502 Noma – 100% Post-consumer recycled polyester fr (fr inside the fiber)

Gabriel tests all polyester fabrics, and tests include all colour options for each fabric. Tests are conducted in accordance with BIFMA's and ACT's recommended cleanability guidelines for use of cleaners, sanitisers and disinfectants on fabrics in hospitals and health care settings. The test result for each colour includes an assessment of the risk for colour change, when bleach is applied to the fabric in the concentrations required in health care environments.

When choosing a bleach-cleanable product, it is important to be aware that a variety of test methods to evaluate bleach resistance exist. Consequently, we recommend that you always ensure that the test method applied to a specific fabric meets the requirements - in terms of bleach concentration, application and contact time - for the specific context and environment in which the fabric will be used.

The test method applied by Gabriel is extremely thorough, and we consider it to be the best test available to assess and inform about the risk for colour change when using chlorine products.

#### Test description

1 ml of hospital grade disinfectant cleaner - diluted in accordance with the manufacturer's instructions - is applied to the centre of the test specimen. The solution is allowed to set for a period of two hours, after which any remaining liquids are blotted up (on both face and back).

The process is repeated for a total of ten times. Two hours after the 10<sup>th</sup> application, three ml of water are applied, excess fluids are blotted up with a clean white cloth, and the test specimen is allowed to air dry. The last step is repeated if chemical residue remains.

The material is evaluated by comparing the test specimen with AATCC Grey Scale for Color change.

#### Rating system – Grades according to AATCC Grey scale

Grade 5 – Very good-excellent Grade 4 – Good Grade 3 – Fair-moderate

Grade 2 – Poor behaviour

Grade 1 – Very poor

Glude I Very poor

#### Acceptance criteria according ACT/BIFMA.

Colour Change:Grade 4 minimumColour Transfer:Not permittedPhysical damage:Not permitted

# Gabriel

Fabric	Colour	Name	Risk for colour changes*	Result
Noma	62155	Light yellow	Low	4
Noma	66283	Light Blue	Low	4
Noma	66285	Dark Blue	Low	4
Noma	61316	Light Brown	Medium	3-4
Noma	68263	Green	Medium	3-4
Noma	67140	Light Blue-Green	Medium	3-4
Noma	66282	Blue	Medium	3-4
Noma	66138	Dark Blue	Medium	3-4
Noma	63154	Light Red Orange	Medium	3-4
Noma	63151	Light Orange	High	3
Noma	63152	Light Orange	High	3
Noma	64286	Orange-Red	High	3
Noma	68324	Medium Green	High	3
Noma	64285	Dark Red	High	1-2

\*) Low risk = Grade 4-5; Medium risk = Grade 3-4; High risk = Grade 3 and below

Gabriel A/S confirms that the above results were obtained after testing the specimen in accordance with the procedures and equipment specified above.

Gabriel A/S

hur legenaard

Kurt Nedergaard Director of CSR & Quality